

Gerd-OFF

Gerd-OFF represents an original approach to GERD management: its components provides a coating to the esophageal mucosa protecting it from both acid (chloridric acid) and non-acid reflux (Pepsine and bile acid) and indirectly stimulating the healing process of the mucosa. Gerd-OFF offers hydration to the upper respiratory tract reducing the atypical symptoms of GERD .

Composition:

- Chondroitin sulfate (CS)
- Hyaluronic acid (HA)
- Bioadhesive carrier poloxamer

Mechanism of action:

- CS provides Anti inflammatory and mucosal protective activity
- HA activates fibroblast ,macrophages and tissue repair
- With CS and HA forms a macromolecular complex with a barrier activity which is prolonged with the esophageal mucosa



Clinical data: 20 patients with non-erosive gastroesophageal reflux symptoms were enrolled in the 2week study by a wash out period of 1-week.

Parameters: Symptom Intensity Difference (SID) , Heartburn and acid regurgitation, Onset of action (Chevrel criteria) and duration of the effect, Tolerability was based on AEs reporting, patient acceptance based on taste and swallowing difficulties possibly linked to syrup viscosity.

Results: Treatment with Gerd-OFF produced a fast relief of GERD symptoms. Secondary analysis showed that significantly more patients receiving Gerd-OFF achieved rapid symptoms disappearance and prolonged symptoms free period.

Indication: Gastro-oesophageal reflux: pyrosis, epigastric pain, acid reflux, irritative cough and dysphonia.

USP: The trio concept of Chondroitin sulfate (CS),Hyaluronic acid (HA) and Bioadhesive carrier poloxamer has proven to provide complete relief from the GERD symptoms

Dose: 10 ml of solution after main meals and at bedtime, or as instructed by your doctor.

Packaging: 10 ml stick packs or bottle (up to 200 ml)

Category: Nutraceutical